



REHASENSE



# PLATFORM

BIG WALKER





# PLATFORM

The Platform is narrow enough to pass through a 70 cm door opening, making it the perfect choice when working in confined spaces and to perform in every situation.

Modern design, rigid and stable but still with a light construction. It has height adjustable gas spring activated arm support.




Clamps for adjusting length & angle of hand grips, arm pad width and locking mechanism for height adjustment. Easy to adjust without tools.




Soft and comfortable PU-pads with optimal thickness for a more stable feeling, Easy width adjustment, Waterproof without seams for better cleaning (hygiene). With lateral support to prevent the arms from sliding off.



Braking system with internal cables underneath. Length and angle adjustable to the users need.



Two parted frame easy to disassemble into a small package for transportation of the device.



The height adjustment lever works in both directions (up and down).



100 mm castors to reach under low objects; all 4 castors directional locking.



A narrow frame allowing it to pass through a 70 cm door opening.

## TECHNICAL DETAILS

Standard colour:

White

Safe weight limit:

150 kg

Arm pads width:

66-80 cm

Arm pads height:

92-132 cm

Length:

80 cm

Height:

18,5 kg

Width between arm pads:

31 - 48 cm

Transportation size

(disassembled in box):

Length x Width x Height:

85 x 77 x 23 cm

Frame material:

Aluminium and steel

Article item:

PLAT1000GB

Without:

- Gas spring support:

PLAT 1000MB

- Brake:

PLAT1000G

- Gas spring support & brake:

PLAT1000M

## ACCESSORIES



One-hand brake



Infusion standard with holder



Oxygen bottle holder

Platform is CE-marked and fulfills the requirements of EU directive 93/42 EEC for medical devices.

Copyright © Rehasense® 2018

Rehasense Benelux  
Markeweg 81  
8398 GM Blesdijke  
Netherlands  
+31 561 450599

Rehasense sp. z o.o  
Sulejowska 45 g  
97-300 Piotrków Trybunalski  
Poland  
+48 516 185 944