



REHASENSE



CLOUD I

ELECTRIC COMMUNITY BED

The **Cloud I** electric community bed, is easy to use, very competitively priced, and its space saving storage facility is of particular interest to Nursing, Care Home and Community environments. The clearly marked handset is simple to use for both carers and end users alike. Despite its unique foldaway design, this bed may be supplied with bed rails for added security.



CLOUD I

ELECTRIC COMMUNITY BED



REHASENSE



TECHNICAL INFORMATION

Product Code:

Max. User Weight Limit:
28 stone / 180 kg

Bed Dimensions:
202 x 102 cm

Overall length with wooden
Head/Foot Boards:
202 cm

Total weight including side
rails and transport kit:
76 kg

Mattress Platform Dimensions:
197 x 88 cm

Max. Height Position:
65.5 cm

Lowest Height Position:
24.5 cm

Bed Weight:
Heaviest of four sections 20 kg

Warranty:
3 years on metal work and
3 years electrical actuator

KEY FEATURES

- Easy assembly and installation
- 75cm twin wheel lockable castors
- Four section bed fully electrically operated
- Electric backrest raises to an angle of 72°
- Electric Height adjustment
- Electric knee break
- Electric Trendelenburg and reverse Trendelenburg functions
- Linak electronics to IPx4 Ingress Protection
- Linak individual function lock-out handset
- Easy clean, low maintenance



CLOUD I BED AND RAIL OPTIONS

- **Cloud I** Bed without Rails
- **Cloud I** Bed with Standard Side Rails 2x2-bar steel full length
- 2x2-bar steel full length side rails (for Community use)
- 3-bar Aluminium wood effect full length side rails
- 3/4 side rails
- Battery back-up option
- wooden bed end panel covers are easily fitted to provide a more "homely" look

EXTRAS

- Bed Lever
- Transport Kit
- Lifting pole
- Bedside cabinet

The **Cloud I** is CE-marked and fulfils the requirements of EU Regulation 2017/745 (MDR) for medical devices.

Copyright © Rehasense® 2021

Rehasense sp. z o.o.
Sulejowska 45g
97-300 Piotrków Trybunalski
Poland