



REHASENSE



SERVER XC-W

OUTDOOR ARTHRITIS ROLLATOR

The **Server XC-W** has the Server XC frame base with the triangle profile which enable a rigid frame and responsive steering. The fore arm supports provide more weight bearing possibilities for the driver who needs more support during walking.



more than mobility

SERVER XC-W

OUTDOOR ARTHRITIS ROLLATOR



REHASENSE


With the wide Server XC base frame the rollator is tested and approved for outdoor use.

It is one of the few Walkers that is folding and can go into a car easily, for transport. With all aluminium construction it is one of the lightest rollators on the market but still supports up to 200 kg user weight. As a side folding design it offers a better walking posture for the user.

The locking carry handle makes sure it stays firmly when folded thus making it easy to transport and store.



Height adjustable forearm supports. Handle with brake lever integrated in the forearm supports adjustable in length.



The handgrips with the brake levers are rotatable to adjust for the best position for the hand.



Hand brakes with parking brake. Very smooth running and can be used with very little force.

TECHNICAL INFORMATION

Colour: Champagne

Max. Users weight: 200 kg

Article no.: NACHL600

Measurements:
(length x width x height)
Open : 77 x 71 x (87 -110) cm
Folded: 77 x 31 x 99 cm

Recommended user height:
140 - 185 cm

Seat height: 62 cm

Seat width: 53 cm

Weight: 12 kg (+ basket 400 g)

Shopping bag:
Width x depth x height:
40 x 17 x 24 cm

Max. load of the bag: 5 kg

Frame material: aluminium

Wheels:
Front: 285 x 44 mm
Rear : 235 x 44 mm
with black soft PU tyres



Stays firmly closed with the lock and is easy to open with the push button.



Easy-to-remove carrying bag with special designed easy-to-go clips.



With a practical small net bag to put between the hand grips for carrying small items.



Large Comfort wheels with black soft PU tyre.



Rotating forearm supports. Can be set in an angle of 22,5° and 45° for maximizing personal comfort.



Safety & relaxation thanks to the Height Adjustable back belt.

Server XC-W is CE-marked and fulfils the requirements of EU Regulation 2017/745 (MDR) for medical devices.

Copyright © Rehasense® 2021

Rehasense sp. z o.o.
Sulejowska 45g
97-300 Piotrków Trybunalski
Poland

Rehasense Xiamen Co., Ltd.
No.289, Wengjiao Road,
Xinyang Industrial Area,
Haicang, Xiamen, China.